Endodontic Micromotor

ENDO-MATE **AT**

ENDO-MATE AT

OPERATION
MANUAL

Thank you for purchasing Endo-Mate AT.

Please read this Operation Manual carefully before use for operating instructions and care and maintenance guidelines so that you can use it for many years to come through a correct use.

Keep this Operation Manual within easy reach of users for future reference.

– Contents –

1. User and Intended Use	1
2. Precautions for handling and operation	1
3. Package Contents	3
4. Part Names	4
5. How to Connect the Parts	6
6. Check before treatment	7
7. Operating Method	7
8. Various Functions	8
9. Maintenance	11
10. Periodical Maintenance Checks	12
11. Error Codes	13
12. Troubleshooting	13
13. Specifications	14
14. Classification of Equipment	15
15. Operation Principle	15
16. Symbol	15
17. Warranty	16
18. Contra Angle Head Correspondence Table	16
19. Spare Parts List	17
20. Disposing product	17
21. EMC Information (Electromagnetic Compatibility Information)	17

1 User and Intended Use

User: Qualified Professionals

Intended Use : Dental Treatment (Endodontics etc.)

2 Precautions for handling and operation

- Please read these precautions carefully and use only as intended or instructed.
- Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.

Class	Degree of Risk
 WARNING	Hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed.
 CAUTION	Hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed.
NOTICE	General product specification information highlighted to avoid product malfunction and performance reduction.

 **WARNING**

- Do not handle the AC adaptor with wet hands. Wet hand contact with electricity may result in an electric shock.
- Keep away from explosive substances and flammable materials.
- If the product overheats or smells of burning, immediately turn off the power and disconnect the AC adaptor. Contact your Authorized NSK Dealer.
- Do not operate close to patients with cardiac pacemakers as there is a danger that it may affect the pacemaker.
- Do not use on the following patients.
 - Those with medical complications or allergies
 - Those who have pre existing conditions (Eg Cardiac, Pulmonary, Renal disturbance or High blood pressure)
 - Those who are pregnant or lactating
 - Patients with cardiac pacemakers and infants
- Be careful not to get water or liquid disinfectant on the Control Unit. This could cause short circuits and lead to fire and/or electric shock.
- The AC Adaptor is the means to cut off commercial power supply. Make sure that the AC Adaptor can be pulled out from the power outlet without delay in an emergency. Do not place any articles within 10cm of the AC Adaptor.

 **CAUTION**


- Read this Operation Manual before use to fully understand the product functions and file for future reference.
- When operating the product always consider the safety of the patient.
- The end user shall be responsible for any judgment that relates to the application of this product to a patient.
- This product does not consider patient's age (except infants), gender, weight or nationality.
- This product does not consider operator's age (mature person), height, weight, gender, or nationality.
- Users are responsible for the operational control, maintenance and continual inspection of this product.
- This device is for indoor use only.
- Keep the Control Unit on a level surface.
- Do not attempt to disassemble the product nor tamper with the mechanism except as recommend by NSK in this Operation Manual.
- Do not allow any impact on to the product. Do not drop the product.
- Operators and all others in the area must wear eye protection and a mask when operating this product.
- Should the product function abnormally, cease operation immediately and contact your Authorized NSK Dealer.
- Do not use high acid water or sterilizing solutions to wipe, immerse or clean the product.
- The Contra Angle Head is delivered in a non-sterile condition and must be autoclaved prior to use.
- Perform regular function and maintenance checks.
- If the product is not used for a long period check it is functioning correctly before using on a patient.
- To avoid clinical downtime it is recommended that a spare be kept on hand in case of a breakdown during surgery.
- Use only the genuine AC Adaptor with this product. Using other adaptors may cause malfunction.
- This product is rated Medical Electrical equipment. EMC (Electromagnetic compatibility) is described in the documentation included.
- Installation and use of this product requires special precautions regarding EMC according to the EMC information.
- Portable and mobile RF communications equipment can affect Medical Electrical equipment. Do not use RF equipment near the product.
- The use of ACCESSORIES such as handpieces and cables, other than those specified by the manufacturer, with the exception of handpieces and cables sold by the manufacturer of this product, as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of this product.

- This product should not be used adjacent to, or stacked with, other equipment. If adjacent or stacked use is necessary, this product should be observed to verify normal operation in the configuration in which it will be used.
- The system may present a possibility of malfunction when used in the presence of an electromagnetic interference wave. Do not install the system in the vicinity of any device which emits magnetic waves. Turn off the Main Power Switch of the system as an ultrasonic oscillation device or an electrode knife is located close to the vicinity of use.
- U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

NOTICE

- No special training is required for this device.
- The micromotor and motor cord during operation could interfere the computers, LAN cables in vicinity of use or could cause noise in radio receivers nearby.


3 Package Contents

No.	Part name	Quantity
1	Control Unit, Motor Handpiece	1
2	AC Adaptor	1
3	Contra Angle Head (iFX75)	1
4	F-type Spray Nozzle	1
5	Handpiece Stand	1

- NOTICE** • The Motor Handpiece cannot be removed from the Control Unit.

4 Part Names

4-1 Keys on the Operation Panel

(1) POWER key

Turns ON/OFF the power for the Control Unit.

(2) SPEED Key (- +)

Sets the rotation speed.

Press (+) or (-) to adjust the speed.

While the key is held down, the value changes continuously. When the value has reached the upper limit or the lower limit, an alarm sounds.

(3) TORQUE Key (- +)


Sets the upper limit of the torque.

Press (+) or (-) to adjust the torque.

While the key is held down, the value changes continuously. When the value has reached the upper limit or the lower limit, an alarm sounds.

- (4) PROGRAM Key (- +)
Press to select the Program No. from 1 to 9.
Press (+) or (-) to change the Program No.
While the key is held down, the Program No. changes continuously.
- (5) REVERSE Key
Press to select the rotation mode (forward rotation, reverse rotation).
- (6) GEAR RATIO Key
Press to select the gear ratio appropriate for the Contra Angle Head in use. (1:1, 4:1, 10:1, 16:1, 20:1)
- (7) AUTO REVERSE Key
Press to select the auto reverse mode.
For details, refer to "8-3 AUTO REVERSE Function".
- (8) MEMORY Key
Press to memorize the values displayed in the LCD Panel in programs.
For details, refer to "8-1 Memory Function".

4-2 LCD Panel

*The diagram shows the Panel when all the lights are on.

- (A) SPEED
Displays the selected rotation speed.
Depending on the speed, **X10** or **X100** may also appear.
- (B) TORQUE
Displays the selected torque upper limit.
- (C) TORQUE BAR GRAPH
The bar graph shows the level of load applied while the Motor Handpiece is rotating.

(D) PROGRAM No.

Displays the selected program No.

(E) GEAR RATIO


Displays the selected gear ratio.

(F) ROTATION MODE

Displays the selected rotation mode. "F" and "R" represent forward rotation and reverse rotation respectively.

(G) AUTO REVERSE MODE

Displays the selected auto reverse mode.

*For details, refer to "8-3 AUTO REVERSE Function".

(H) ALARM SOUND VOLUME

Displays the selected alarm sound volume.

*For details, refer to "8-5 Selecting the Alarm Sound Volume".

5 How to Connect the Parts

5-1 Connection of the AC Adaptor

- 1) Position the AC Adaptor plug so that ➔ mark is on its bottom face, then insert it to the AC Adaptor inlet on the side of the Control Unit. (Fig.1)
- 2) Plug the AC Adaptor into a power outlet for medical use.


Fig. 1

CAUTION

- Do not pull out the AC Adaptor by holding and pulling its cords.
- Make sure that the device has stopped before connecting or disconnecting the AC Adaptor.

5-2 Mounting of the Contra Angle Head

(1) Mounting

- 1) Make sure that the power is OFF.
- 2) Push the Contra Angle Head till its alignment pins click into the alignment slots in the Motor Handpiece. The Contra Angle Head can be mounted at 6 patterns of angles. (Fig.2)

(2) Removing

- 1) Make sure that the power is OFF.
- 2) Pull out the Head straight without twisting.


Fig. 2

**CAUTION**

- When mounting and removing the Contra Angle Head, do not twist or turn.
- Do not mount heads other than the Contra Angle Head manufactured and specified by NSK.
*Refer to "18. Contra Angle Head Correspondence Table".

5-3 Attaching the file**5-3-1 Mounting**

- 1) Insert the file until it is correctly seated in place.
- 2) Depress the Push Button and insert the file into the chuck until the file "notch" mechanism engages. Release the button.
- 3) Ensure that the file is secure by gently pulling and pushing the file WITHOUT depressing the Push Button.


Fig. 3

5-3-2 Removing

Depress the Push Button and remove the file.

**CAUTION**

- Make sure that the power is OFF before attaching or removing the file.
- Do not mount or remove the file until the drive motor has completely stopped.
- Depressing the push button while the handpiece is rotating will result in OVERHEATING of the handpiece head and burn injury to the patient.
- Always keep the file shank clean. Debris inside the chuck could reduce file concentricity.
- DO NOT use bent, deformed, cracked or non-ISO standard files as the file may break or bend during use.
- Do not exceed the drive motor speed or torque recommended by the file manufacturer.

6 Check before treatment

Check that the Head Cap is firmly tightened. Also check for handpiece vibration, noise and overheating. If any abnormalities are found do not use the handpiece and contact your Authorized NSK Dealer.

7 Operating Method**7-1 Preparation**

- 1) Press down the POWER key for one second or longer to turn ON the power.
- 2) Press the PROGRAM Key to select the Program No. to use.
- 3) Press the GEAR RATIO Key to select the gear ratio appropriate for the Contra Angle Head mounted.
- 4) Press the SPEED Key to set the speed, the REVERSE Key to set the rotation mode, and the TORQUE Key to set the torque upper limit.
- 5) Press the AUTO REVERSE Key to select the Auto Reverse mode. (Refer to "8-3 AUTO REVERSE Function")

7-2 Operation

7-2-1 Alternate operation

- 1) When the ON/OFF Key is pressed briefly, the Motor Handpiece starts rotation.
- 2) When the ON/OFF Key is pressed again, the rotation stops.

7-2-2 Intermittent operation

- 1) Press and hold down the ON/OFF Key for one second or longer.
- 2) While the Key is held down, the Motor Handpiece rotates.
- 3) When the ON/OFF key is released, the rotation stops.


Fig. 4

For fine-tuning the rotation speed and the torque upper limit, press the SPEED Key and the TORQUE Key to adjust the values.

When the treatment is over, press and hold down the POWER Key for one second or longer to turn off power.

8 Various Functions

8-1 Memory Function

The values (speed, torque, gear ratio, auto reverse mode) set as required can be saved under a chosen program.

- 1) Press the PROGRAM Key to select a Program No. for saving the setting into memory.
- 2) Set the speed, torque, gear ratio, auto reverse mode as required.
- 3) Press and hold down the MEMORY Key for one second or longer. When the alarm sounds, the setting has been saved in the selected program.

- NOTICE**
- A program cannot be saved while the Motor Handpiece is rotating.
 - The Rotation Modes cannot be saved in a program. Also, when the rotation mode is set at "R (reverse)", the setting cannot be saved.
 - When the PROGRAM Key is pressed and held down, Program No. changes continuously.
 - The MEMORY Key needs to be pressed and held down for two seconds or longer to save the program. If the MEMORY Key is not pressed long enough, the setting is not saved in the memory. If the Program No. is changed using the PROGRAM Key before the MEMORY Key is pressed, the setting remains the original saved values (cancel function).

8-2 Restoring the Default Setting of Programs

The factory setting can be restored as follows.

- 1) Turn OFF the power.
- 2) While pressing the AUTO REVERSE Key and the MEMORY Key together, press and hold down the POWER key for one second or longer.
- 3) When the alarm sounds and "5 8 8" has appeared on the LCD Panel, release both keys. Then while "5 8 8" is displayed, press the MEMORY Key.
- 4) When "F 17" appears on the LCD Panel, the default setting has been restored.

- NOTICE** • This function deletes all the settings stored in the programs and restores the original values (factory default settings). If necessary, note down the current settings before using this function.

8-3 AUTO REVERSE Function

While the Motor Handpiece is rotating, the alarm sounds when the torque has reached half of the set upper level. (At this point, the bar graph looks like this. )

When the torque gets close to the specified upper limit, the alarm sound changes. (At this point, the bar graph looks like this. )

If load is applied further and the torque exceeds the set upper limit, the Motor Handpiece starts operating in the selected Auto Reverse mode. The Auto Reverse mode can be selected from among the three modes shown below.

AUTO REVERSE

When the torque reaches the set load level, the forward rotation of the Motor Handpiece stops, and the reverse rotation starts. When load is no longer applied, the forward rotation resumes automatically.


Fig. 5

AUTO STOP

When the torque reaches the set load level, the forward rotation of the Motor Handpiece stops, and the reverse rotation starts. When load is no longer applied, the rotation stops. To resume the rotation (forward rotation), press the ON/OFF Key.


Fig. 6

AUTO REVERSE OFF

When the torque has reached the set load level, the rotation stops. To resume the rotation (forward rotation), press the ON/OFF Key.

- NOTICE**
- When "R (reverse)" is selected as the ROTATION MODE, the AUTO REVERSE function does not work.
 - When load is applied continuously for a long time, the rotation may stop automatically to prevent overheating. (In such a case, "EH" is displayed in the TORQUE on the LCD Panel.) If this happens, stop using the Motor Handpiece and wait till it cools down.

8-4 Calibration

This function reduces variance in the rotation speed among individual Motor Handpieces and variance in the torque among individual Contra Angle Heads.

- 1) Remove the file.
- 2) Lubricate the Contra Angle Head.
- 3) Turn ON the power.
- 4) Press both (-) and (+) of the SPEED Key together for one second or longer.
- 5) The alarm sounds and "EH" appears on the LCD Panel.
- 6) Mount the Contra Angle Head onto the Motor Handpiece.
- 7) Press the MEMORY Key.
- 8) The rotation of the Motor Handpiece starts. Leave it until the rotation stops.
- 9) The rotation stops and the alarm sounds. When the display on the LCD Panel returns to the original state, the calibration is complete.
- 10) To stop the calibration, turn OFF the power.

- NOTICE**
- Make sure to lubricate the Contra Angle Head prior to calibration. If the Head is not lubricated, the calibration cannot be performed properly.
 - Do not touch the rotational axis of the Contra Angle Head or apply load during calibration. The measuring may not be accurate.
 - This function cannot absorb individual differences of Motor Handpieces and Contra Angle Heads completely.

8-5 Selecting the Alarm Sound Volume

The alarm sound volume can be selected from High, Low and Limited OFF.

- 1) Press and hold down both (-) and (+) of the PROGRAM Key together.
- 2) The ALARM SOUND VOLUME display on the LCD Panel changes to represent the sound volume.
- 3) When the volume of your choice is displayed, release the PROGRAM Key to set the sound volume.

High

Low

Limited OFF

*The alarm sounds in the low volume for confirmation and error alert. The alarm does not sound during reverse rotation or when the torque has reached the specified upper limit.

- NOTICE**
- The alarm setting is saved even when the power is switched OFF.

8-6 Auto Power Off Function

After 10 minutes of non-activity in key operation, power is switched OFF automatically for saving energy and preventing an accidental operation. When the Motor Handpiece is rotating, however, the power is not switched off even if keys are not operated.

9 Maintenance

After each patient maintain the product as follows.

9-1 Cleaning (Contra Angle Head)

- 1) Remove debris from the product. Do not use a wire brush.
- 2) Wipe clean with alcohol-immersed cotton swab or cloth.

 This icon denotes that the product can be washed via Thermo Disinfector. (IFX75)
Refer to the Thermo-Disinfector manual.

CAUTION

- After washing with Thermo-Disinfector and prior to lubrication, dry the product until all internal moisture is thoroughly removed. Thermo-Disinfector moisture remaining inside the product could reduce the effect of lubrication and could cause corrosion inside of the product. (IFX75)
- To clean the product never use any solvent such as benzene or thinner.

9-2 Cleaning (Motor Handpiece, Control Unit)

- 1) Turn OFF the power.
- 2) Detach the AC Adaptor.
- 3) Wipe clean the surface, first with a damp cloth, then an alcohol-soaked cloth.

CAUTION

- To clean the product never use any solvent such as benzene or thinner.

9-3 Lubrication (Contra Angle Head)

■ NSK PANA SPRAY Plus

Apply NSK PANA SPRAY Plus every time after each use and/or before autoclaving.

- 1) Remove the file from the Contra Angle Head.
- 2) Screw the F-Type Spray Nozzle into the nozzle portion of PANA SPRAY Plus by turning it 10 times or so.
- 3) Insert the spray nozzle into the rear of the Contra Angle Head. Hold the Contra Angle Head and spray for 2 to 3 seconds till oil comes out of the tip of the Contra Angle Head. Repeat lubrication until the oil coming out of the tip is clear of debris, etc. (Fig. 7)


Fig. 7

CAUTION

- Only the Contra Angle Head may be lubricated.
- When applying spray be sure to firmly hold the Contra Angle Head to prevent the Contra Angle Head from slipping out of the hand due to the spray pressure.
- Hold the spray can upright.

9-4 Sterilization (Contra Angle Head)

Sterilize the Contra Angle Head by autoclave sterilization. Remove the file after each patient and sterilize as noted below.

1) Insert into an autoclave pouch. Seal the pouch.

2) Autoclavable under the conditions below.

Autoclave for more than 20 min. at 121°C, or 15 min. at 132°C, or 3 min. at 134°C.

3) The Contra Angle Head should remain in the autoclave pouch until required for use.

CAUTION

- Only the Contra Angle Head may be sterilized.
- Do not autoclave the product with other instruments even when it is in a pouch. This is to prevent possible discoloration and damage to the product from chemical residue on other instruments.
- Keep the product in suitable atmospheric pressure, temperature, humidity, ventilation, and sunlight. The air should be free from dust, salt and sulphur.
- Immediately after use, the product should be cleaned, lubricated and sterilized. If blood remains on the external or internal surfaces it can become clotted and cause rust.
- Do not heat or cool the product too quickly. Rapid change in temperature could cause damage to the product.
- If the sterilizer chamber temperature may exceed 135°C during the dry cycle then delete the dry cycle.
- Autoclave sterilization is recommended for the product. The validity of other sterilization methods is not confirmed.
- Do not touch the product immediately after autoclaving as it will be very hot and must remain in a sterile condition.

NOTICE • NSK recommends Class B sterilizers as stated in EN13060.

10 Periodical Maintenance Checks

Perform periodical maintenance checks every three months, referring to the check sheet below. If any abnormalities are found, contact your Authorized NSK Dealer.

Points to check	Details
Looseness of the Head Cap	Check that the Head Cap is not loose.
Rotation	Rotate the handpiece and check for abnormalities such as abnormal rotation, vibration, noise, and overheating.
Display	Immediately after the power is turned On, all the displays lights on the LCD Panel should be on. Check if any of the displays lights are missing.

11 Error Codes

When the device has stopped due to failure, excessive load, disconnection and incorrect usage, the error code is displayed in the Display Panel of the Control Unit. When this happens, press the ON/OFF Key to repeat the error check. If no problems are detected, the error is cancelled and operation can be resumed. If the error is displayed again, refer to the table below and take appropriate actions.

<When the Motor Handpiece is rotating>

Error code	Error	Cause	Remedy
$\mathcal{E} - 0$	Self-check error	Circuit failure	Contact your NSK Dealer.
$\mathcal{E} - 1$	Excess current	Motor Handpiece is locked (in the Auto Reverse mode)	Stop applying load.
$\mathcal{E} - 2$	Excess voltage	Circuit failure	Contact your NSK Dealer.
$\mathcal{E} - 4$	Motor overheating	A high load has been continuously applied on the Motor Handpiece for a relatively long time.	Let the motor cool down before resuming use.
$\mathcal{E} - 5$	Brake circuit error	Circuit failure	Contact your NSK Dealer.
$\mathcal{E} - 6$	Locked axis	The Motor Handpiece is locked (in the start-up)	Stop applying load.

<In calibration>

Error code	Error	Cause	Remedy
$\mathcal{E} \mathcal{E} 0$	Above the upper limit	Motor Handpiece or Contra Angle Head has come to end of its life.	Replace the Contra Angle Head. If the error persists, the Motor Handpiece may have come to the end of its life. Contact your Dealer.
$\mathcal{E} \mathcal{E} 1$	Below the lower limit		

12 Troubleshooting

When a problem is detected, check the following again before requesting a repair.

If none of these is applicable or if the trouble is not remedied even after an action has been taken, a failure of this product is suspected. Contact your Authorized NSK Dealer.

<Control Unit>

Trouble	Cause	Remedy
The power does not turn on.	AC Adaptor is not connected to the Control Unit.	Check the connection.
	The plug of AC Adaptor is not connected to the power outlet. Or, electricity is not supplied to the power outlet.	Check the connection.
	The internal fuse has been activated.	Contact your NSK Dealer.

<Motor Handpiece>

Trouble	Cause	Remedy
Motor Handpiece does not rotate.	Open circuit inside the Motor Handpiece or the cord of Motor Handpiece.	Contact your NSK Dealer.
Motor Handpiece does not rotate. (Error code " $\mathcal{E} - 1$ " is displayed.)	Contra Angle Head is clogged up.	Clean or replace the Contra Angle Head.
	Short circuit inside the Motor Handpiece or the cord of Motor Handpiece.	Contact your NSK Dealer.

The Motor Handpiece does not rotate. (" - - -" and the speed are displayed alternately.)	Contra Angle Head is clogged up.	Clean or replace the Contra Angle Head.
	Short circuit inside the Motor Handpiece or the cord of Motor Handpiece.	Contact your NSK Dealer.
When the power is turned On, the alarm sounds and the Motor Handpiece does not rotate.	The ON/OFF button is pressed when the power is turned ON.	Check the ON/OFF button.
	Short circuit inside the ON/OFF button.	Contact your NSK Dealer.

13 Specifications

<Control Unit, Motor Handpiece>

Model	NE305
Rated input	DC20V 0.5A
Rated output	DC7V 0.4A
Dimensions (Control Unit)	W268 x D230 x H103mm
Dimensions (Motor Handpiece)	ø20×L108.5mm
Weight	342g (including the Motor Handpiece)

<AC Adaptor>

Model	NE169	NE169-02	NE169-03
Rated input	AC100V 50/60Hz	AC120V 60Hz	AC230V 50/60Hz
Cord length	Plug side 1.8m/Unit side 1.8m		Plug side 2m/ Unit side 1.8m
Dimensions	W70 x D110 x H58mm		
Weight	918g	933g	974g

<Contra Angle Head>

Model	iFX75
Max. Rotation Speed (Motor)	13,000min ⁻¹
Max. Rotation Speed (Handpiece)	813min ⁻¹
Gear Ratio	16:1 Reduction
File Type	ISO 1797-1 Type1 ø2.35mm Engine File
Chuck Length	9.4mm
File Motion	360°Rotation
Weight	41g

	Temperature	Humidity	Atmospheric Pressure
Use Environment	10 - 40°C	10 – 85%	-
Transportation and Store Environment	-10 - 50°C	10 – 85%	500 - 1060hPa

* No Condensation

* Using in conditions outside the limits above may cause malfunction.

14 Classification of Equipment

- Type of protection against electric shock :
 - Class II equipment 
- Degree of protection against electric shock :
 - Type B applied part  (Applied Part: File)
- Method of sterilization or disinfection recommended by the manufacturer :
 - See "9-4 Sterilization (Contra Angle Head)"
- Degree of safety of application in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide :
 - Equipment NOT suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- Mode of operation :
 - Continuous operation

15 Operation Principle

This system consists of the Contra Angle Head, the Control Unit with the Motor Handpiece, and the AC Adaptor. By operating the ON/OFF Key of the Motor Handpiece, power is supplied to the motor and the file attached to the Contra Angle Head rotates.

16 Symbol

This product is Autoclavable up to Max. 135°C.

This product can be washed via Thermo Disinfectant.

Conforms to CE European Directive of "Medical device directive 93/42/EEC."

Manufacturer.

Authorized representative in the European community.

Follow the waste of electric and electronic equipment (WEEE) Directive (2002/96/EC) for product and accessory disposal.

Consult operation instructions.

Caution, Refer to attached instructions

Class II equipment.

Type B applied part.

Short circuit proof safety isolating transformer.

For indoor use only.

To identify fuse boxes or their location.

This product meets UL safety standard requirements.

This mark is machine tool, electrical equipment and medical equipment Safety mark in Canada.

Product safety electrical and material.

Marking on the outside of Equipment or Equipment parts that include RF transmitters or that apply RF electromagnetic energy for diagnosis or treatment.

TUV Rhineland of North America is a Nationally Recognized Testing Laboratory (NRTL) in the United States and is accredited by the Standards Council of Canada to certify electro-medical products with Canadian National Standards.

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

17 Warranty

NSK products are warranted against manufacturing errors and defects in materials. NSK reserves the right to analyze and determine the cause of any problem. Warranty is voided should the product be not used correctly or for the intended purpose or has been tampered with by unqualified personnel or has had non NSK parts installed. Replacement parts are available for seven years beyond discontinuation of the model.

18 Contra Angle Head Correspondence Table

Model	Order Code	Gear Ratio	
iFX75	C1098	16:1 Reduction	-Mini head -360° rotation -For the Ni-Ti file (ø2.35)
MP-F20R	Y110044	20:1 Reduction	
MP-F4R	Y110047	4:1 Reduction	
NML-F20R	Y110053	20:1 Reduction	
NML-F16R	Y110019	16:1 Reduction	
NML-F4R	Y110018	4:1 Reduction	
MPA-F20R	Y110045	20:1 Reduction	-Mini head -360° rotation -For Ni-Ti file (ø2.35) -For connection to Apex Locator
MPA-F16R	Y110097	16:1 Reduction	
MPA-F4R	Y110051	4:1 Reduction	
MPAS-F20R	Y110046	20:1 Reduction	
MPAS-F16R	Y110098	16:1 Reduction	
MPAS-F4R	Y110052	4:1 Reduction	

19 Spare Parts List

Model	Order Code
F-type Spray Nozzle	Z256090

20 Disposing product

In order to avoid the health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose the product for you.

21 EMC Information (Electromagnetic Compatibility Information)

Guidance and manufacturer's declaration - Electromagnetic Emissions			
The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that is used in such an environment.			
Emissions test	Compliance	Electromagnetic environment - guidance	
RF emissions CISPR11/EN55011	Group 1	The product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF emissions CISPR11/EN55011	Class B	The product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonic emissions EN/IEC61000-3-2	Class A		
Voltage fluctuations/flicker emissions EN/IEC 61000-3-3	Applicable		
Guidance and manufacturer's declaration - Electromagnetic Immunity			
The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that is used in such an environment.			
Immunity test	EN/IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) EN/IEC61000-4-2	±(2, 4) 6kV contact ±(2, 4) 8kV air	±(2, 4) 6kV contact ±(2, 4) 8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst EN/IEC61000-4-4	±2kV for power supply lines ±1kV for input/output lines	±2kV for power supply lines ±1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN/IEC61000-4-5	±1kV line(s) to line(s) ±2kV line(s) to earth	±1kV line(s) to line(s) ±2kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines EN/IEC61000-4-11	<5% Ut (>95% dip in Ut) for 0.5 cycles 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 0.5 cycles <5% Ut (<95% dip in Ut) for 5 sec	<5% Ut (>95% dip in Ut) for 0.5 cycles 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 0.5 cycles <5% Ut (<95% dip in Ut) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the product requires continued operation during power mains interruptions, it is recommended that the product be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz)magnetic field EN/IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: 'Ut' is the AC mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration - Electromagnetic Immunity			
The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that is used in such an environment.			
Immunity test	EN/IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF EN/IEC61000-4-6	3Vrms 150kHz to 80MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the product, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80MHz to 800MHz $d = 2.3\sqrt{P}$ 800MHz to 2.5GHz Where d is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, and (d) is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters as determined by an electromagnetic site survey ⁽⁶⁾ should be less than the compliance level in each frequency range ⁽⁶⁾ . Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF EN/IEC61000-4-3	3V/m 80MHz to 2.5GHz	3 V/m	
NOTE 1 At 80MHz and 800MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a: Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the product is used exceeds the applicable RF compliance level stated above, the product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the product.			
b: Over the 150kHz to 80MHz frequency range, the field strength should be less than 3V/m.			

Cables and accessories	Maximum length	Complies with	
Handpiece cords	1.5m (Unshielded)	RF emissions, CISPR11	Class B/ Group 1
Power cords Plug side	1.8m/2.0m (Unshielded)	Electrostatic discharge (ESD):	EN/IEC61000-4-2
Control Unit side	1.8m (Unshielded)	Electrical fast transient/burst:	EN/IEC61000-4-4
		Surge:	EN/IEC61000-4-5
		Voltage dips, short interruptions and voltage variations on power supply input lines:	EN/IEC61000-4-11
		Power frequency(50/60Hz) magnetic field:	EN/IEC61000-4-8
		Conducted RF:	EN/IEC 61000-4-6
		Radiated RF:	EN/IEC61000-4-3

Recommended separation distances between portable and mobile RF communications equipment and the product			
The product is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the product can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the product as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d = 1.2\sqrt{P}$	80MHz to 800MHz $d = 1.2\sqrt{P}$	800MHz to 2.5GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance 'd' in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where 'P' is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE1: At 80MHz and 800MHz, the higher frequency range applies.			
NOTE2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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NSK United Kingdom Ltd
NSK America Corp
NSK Oceania Pty Ltd
NSK Asia Pte Ltd